

LCG - BEII 2016 - Introduction to the ChIP-seq technology

Jacques van Helden

2016-02-10

Schema of a ChIP-seq experiment

1. cross-link DNA with protein (covalent linkage with formaldehyde)
2. DNA fragmentation by sonication
3. immunoprecipitation
4. unlink DNA from protein
5. fragment size selection (typically ~300bp)
6. sequencing

Typical result

- ▶ A hard drive with huge files
 - ▶ one or several files per sample (one per run)
 - ▶ several 10s of million reads per file
 - ▶ read length: typically between 50 and 150bp
 - ▶ Sequence format: fastq
 - ▶ Quality index associated to each nucleotide of each read

Example of fastq file

```
@SRR576933.1 HWUSI-EAS1789_0000:2:20:1269:14140/1
AAGCATGGAATAACCGCCTGGTGAATGCTCGCCATA
+
dcd`\dddddaeacecdac`c\cca`bTbbdddYd_
@SRR576933.2 HWUSI-EAS1789_0000:2:20:1270:19579/1
TGGAGGCTGACCACGATAAGCTGCCGCTGGTGGTGC
+
dceYc^\cddd^dddTccc`daYdbdaad`]``XTU
@SRR576933.3 HWUSI-EAS1789_0000:2:20:1270:17351/1
AGTGCGATGCCGTTACCCGGTTTTCTTTATCATTA
+
dddddcc\cc^`c\cddadcdadbbbc]]]aa^ddT
...
```

Analysis steps

- ▶ **read quality**: checking the quality of the raw short reads
- ▶ **read mapping**: aligning the reads against the reference genome, to identify their genomic location
- ▶ **peak-calling**: identification of genomic regions that are enriched in reads in immunoprecipitated samples versus background (genomic input, mock, ...)

Peak-calling comparison workflow





Figure 1: Peak-calling comparison workflow. Various software tools and parameters were tested on the same ChIP-seq dataset (Abd-A TF in *Drosophila melanogaster*). Analysis and figure by Claire Rioualen, 2016.



From read to peaks

References